

List of Chemistry Laboratory Equipments

Apparatus

Apparatus	Qty
1. Glass rod	12
2. Laboratory Thermometer	6
3. Burette (50 ml)	3
4. R.B. flask (250 ml)	15
5. Volumetric flask (500 ml)	12
6. Volumetric flask (100 ml)	2
7. R.B. flask (500 ml)	6
8. Beaker (1000 ml)	6
9. Beaker (500 ml)	11
10. Funnel (75 mm)	11
11. Separating funnel (250mm)	9
12. Conical flask (150 mm)	6
13. Chemical store	28
14. Beaker (250 ml)	34
15. Conical flask (250 ml)	23
16. Chemical store	28
17. Beaker (1000 ml)	23
18. China dish	16
19. Bunsen	8
20. Test tube stand	15
21. Test tube	57
22. Measuring cylinder (500ml)	1
23. Measuring cylinder (250ml)	6
24. Measuring cylinder (100ml)	4
25. Measuring cylinder (50 ml)	1
26. Burette stand	9
27. Tripod stand	8
28. Wire gauge	8
29. Spatula	9
30. Test tube holders	4
31. Iron Tongs	6

List 2: Chemicals

1. Ammonium Carbonate
2. Ammonium Ferrous Sulphate
3. Ammonium Chloride
4. Acetomine
5. Ammonia
6. Aluminium Potassium Sulphate

7. Aluminium Sulphate
8. Acetic Acid
9. Acetone
10. Boric Acid
11. Barium Chloride
12. Benzene
13. Calcium Carbonate
14. Caustic Soda
15. Copper Nitrate
16. Calcium Oxide
17. Ethanol
18. Ethanol
19. Ferrous Sulphate Heptahydrate
20. Ferrous Sulphate Fused Stick
21. Ferrous Sulphate Crystalline
22. Ferrous Sulphate
23. Ferrous Sulphate
24. Ferric Chloride
25. Hydrochloric Acid
26. Hydrochloric Acid
27. Litmus Blue
28. Litmus Red Indicator
29. Lead Acetate
30. Litmus Blue Indicator
31. Lead Acetate
32. Leishman Stain
33. Methanol
34. Methylene Blue
35. Nitric Acid
36. Naphthalene
37. Oxalic Acid
38. Oxalic Acid
39. Potassium Permanganate
40. Potassium Dichromate
41. Potassium Hydroxide
42. Potassium Sulphate
43. Potassium Permanganate
44. Potassium Sulphate
45. Potassium Permanganate
46. Phenolphthalein
47. Potassium Iodide
48. Sodium Hydrogen Carbonate
49. Sodium Hydroxide Nicks
50. Sulphuric Acid
51. Sodium Carbonate
52. Starch Soluble
53. Sodium Nitrate

(Signature)
Principal
Jyoti International School
Vijay Garh Road, Berha Patti,
Hodal-121106 (Hr.)
Affiliation No. 531888